



RAW

- **Hide Unused Fields:** If there is not a specification for a test, then the field will be hidden from view.
- **Automatically Include Calibrations:** When printing Antibiotic Reports, calibrations will automatically be included on the report.

- Raw Menu > Raw Defaults > Screen Color & Settings > **Hide Unused Fields.**
- Raw Menu > Raw Defaults > Test Defaults > **Automatically Include Calibrations.**



FLUID

- **Warn on Code Dates:** If the code differs from what is set-up, an alert box will pop up. This can be turned off also.
- **Outside Lab Comparison Reports:** There are three reports that will compare In-house lab testing with external lab testing results.

- Fluid Menu > Fluid Defaults > Data Entry Defaults > **Warn on Code Dates.**
- Fluid Menu > Fluid Reports > In-house Reports > **Outside Lab Comparison Reports.**



ICE CREAM

- **Four Search By Dates:** There are now four different date ranges to Search By. Entered Date is the date entered into QA Studio which may differ from the Test Date.
- **Packaging Information:** There is a field for Packaging in the main window. The sequence of fields including the Packaging field have been rearranged for faster data entry.

- Ice Cream Menu > Ice Cream Window > Find > **Search By.**
- Ice Cream Menu > Ice Cream Window > Flavoring Data Tab > **Packaging.**



CULTURE

- **Minutes to Keep Initials:** There is an option to keep initials in the main window for a specific length of time.
- **Metal Detectors & Viscosity Tests:** There are set-ups for three metal detectors and three viscosity tests. Each field can be named in the Test Defaults window.

- Culture Menu > Culture Defaults > Culture Defaults > **Minutes to Keep Initials.**
- Culture Menu > Culture Defaults > Test Defaults > **Metal Detectors & Viscosity Tests.**

Workshops offered
by Randolph
Associates, Inc.

Dairy Technology Workshop
September 24 - 27, 2007
Birmingham, AL

HTST Workshop
December 3 - 5, 2007
Murfreesboro, TN



Authorized Dealer of QA Studio

QA Studio is committed to meeting and exceeding your needs for a dairy laboratory software program.

We would like your feedback. Please email us at support@gastudio.com